imprimir

INFORMAÇÕES

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Referências

1. **Maki DG.** Infections caused by intravascular devices used for infusion therapy: pathogenesis, prevention, and management. In: Bistrot AI, Waldvogel FA, eds. *Infections Associated with Indwelling Medical Devices*. 2d ed Washington, DC: ASM Pr; 1994:155-205.

2. **Peters G, Locci R, Pulverer G.** Microbial colonization of prosthetic devices. II. Scanning electron microscopy of naturally infected intravenous catheters. *Zentralbl Bakteriol Mikrobiol Hyg [B]*. 1981; 173:293-9.

3. **Maki DG, Ringer M, Alvarado CJ.** Prospective random sed trial of povidone-iodine, alcohol, and chlorhexidine for prevention of infection associated with central venous and arterial catheters. *Lancet*. 1991; 338:339-43.

4. **Levin A, Mason AJ, Jindal KK, Fong IW, Golstein MB.** Prevention of hemodialysis subclavian vein catheter infections by topical povidone-iodine. *Kidney Int*. 1991; 40:934-8.

5. **Maki DG, Cobb L, Garman JK, Shapiro JM, Ringer M, Helgerson RB.** An attachable silver-impregnated cuff for prevention of infection with central venous catheters: a prospective randomized multi-center trial. *Am J Med*. 1988; 85:307-14.

6. **Lowbury EJ, Lilly HA.** Use of 4 per cent chlorhexidine detergent solution (Hibiscrub) and other methods of skin disinfection. *Br Med J*. 1973; 1:510-5.

7. **Maki DG, Hecht J.** Antiseptic-containing handwashing agents reduce nosocomial

infectors-a prospective study [Abstract]. In: Program and Abstracts of the Twenty-second Interscience Conference on Antimicrobial Agents and Chemotherapy, October 1982, Miami Beach, FL. Washington, DC: American Society for Microbiology; 1982:699.

8. Doebbeling BN, Stanley GL, Sheetz CT, Pfaller MA, Houston AK, Annis L, et al. Comparative efficacy of alternative hand-washing agents in reducing nosocomial infections in intensive care units. *N Engl J Med.* 1992; 327:88-93.

9. Loe H. Chlorhexidine in the prophylaxis of dental diseases. Introduction. *J Periodontal Res.* 1973; 8(Suppl 12):5-6.

10. Crosfill M, Hall R, London D. The use of chlorhexidine antiseptics in contaminated surgical wounds. *Br J Surg.* 1969; 56:906-8.

11. Fabbri L, Grimaldi C, Zucchelli P. Peritonitis in CAPD: treatment with chlorhexidine. *Dialysis Transplant.* 1982; 11:483-6.

12. Ball AJ, Carr TW, Gillespie WA, Kelly M, Simpson RA, Smith PJ. Bladder irrigation with chlorhexidine for the prevention of urinary infection after transurethral operations: a prospective controlled study. *J Urol.* 1987; 138:491-4.

13. Christensen KK, Christensen P, Dykes AK, Kahlmeter G. Chlorhexidine for prevention of neonatal colonization with group B streptococci. III. Effect of vaginal washing with chlorhexidine before rupture of the membranes. *Eur J Obstet Gynecol Reprod Biol.* 1985; 19:231-6.

14. Clarke AM. Topical use of silver sulfadiazine and chlorhexidine in the prevention of infection in thermal injuries. *Med J Aust.* 1975; 1:413-9.

15. Inman RJ, Snelling CF, Roberts FJ, Shaw K, Boyle JC. Prospective comparison of silver sulfadiazine 1 per cent plus chlorhexidine digluconate 0.2 per cent (Silvazine) and silver sulfadiazine 1 per cent (Flamazine) as prophylaxis against burn wound infection. *Burns Incl Therm Inj.* 1984; 11:35-40.

16. Carr HS, Wlodkowski TJ, Rozenkranz HS. Silver sulfadiazine: in vitro antibacterial activity. *Antimicrob Agents Chemother.* 1973; 4:585-90.

17. Hoffmann S. Silver sulfadiazine: an antibacterial agent for topical use in burns. A review of the literature. *Scand J Plast Reconstr Surg.* 1984; 18:119-26.

18. Sawhney CP, Sharma RK, Rao KR, Kaushish R. Long-term experience with 1 per cent topical silver sulphadiazine cream in the management of burn wounds. *Burns.* 1989; 15:403-6.

19. Koo DS, Zhen S, Zhen ZD, Shi XW, Xiang SJ. Assessment of topical therapy of the burn wound with silver sulphadiazine after its use for 15 years in a burn unit. *Burns.* 1989; 15:193-6.

20. Mermel LA, McCormick RD, Springman SR, Maki DG. The pathogenesis and epidemiology of catheter-related infection with pulmonary artery Swan-Ganz catheters: a prospective study utilizing molecular subtyping. *Am J Med.* 1991; 91:197S-205S.

21. Maki DG, Ringer M. Risk factors for infusion-related phlebitis with small peripheral venous catheters. A randomized controlled trial. *Ann Intern Med.* 1991; 114:845-54.

22. Knaus WA, Draper EA, Wagner DP, Zimmerman JE. APACHE II: a severity of disease classification system. *Crit Care Med.* 1985; 13:818-29.

23. **Mermel LA, Stolz SM, Maki DG.** Surface antimicrobial activity of heparin-bonded and antiseptic-impregnated vascular catheters. *J Infect Dis.* 1993; 165:920-4.
24. **Maki DG, Weise CE, Sarafin HW.** A semiquantitative culture method for identifying intravenous-catheter-related infection. *N Engl J Med.* 1977; 296:1305-9.
25. **Lennette E, Balows A, Hausler W, Shadomy H, eds.** Manual of Clinical Microbiology. 3d ed. Washington, DC: American Society for Microbiology; 1980.
26. **Alvarado CJ, Stolz SM, Maki DG.** Nosocomial infections from contaminated endoscopes: a flawed automated endoscope washer. An investigation using molecular epidemiology. *Am J Med.* 1991; 91(Suppl 3B):272S-80S.
27. **Maki DG.** Risk factors for nosocomial infection in intensive care. "Devices vs nature" and goals for the next decade. *Arch Intern Med.* 1989; 149:30-5.
28. **Bozzetti F, Terno G, Bonfanti G, Scarpar D, Scotti A, Ammatuna M, et al.** Prevention and treatment of central venous catheter sepsis by exchange via a guidewire. A prospective controlled trial. *Ann Surg.* 1983; 198:48-52.
29. **Gregory JA, Schiller WR.** Subclavian catheter changes every third day in high risk patients. *Am Surg.* 1985; 51:534-6.
30. **Hagley MT, Martin B, Gast P, Traeger SM.** Infectious and mechanical complications of central venous catheters placed by percutaneous venipuncture and over guidewires. *Crit Care Med.* 1992; 20:1426-30.
31. **Michel LA, Bradpiece HA, Randour P, Pouthier F.** Safety of central venous catheter change over guidewire for suspected catheter-related sepsis. A prospective randomized trial. *Int Surg.* 1988; 73:180-6.
32. **Snyder RH, Archer FJ, Endy T, Allen TW, Condon B, Kaiser J, et al.** Catheter infection. A comparison of two catheter maintenance techniques. *Ann Surg.* 1988; 208:651-3.
33. **Eyer S, Brummit C, Crossley K, Siegel R, Cerra F.** Catheter-related sepsis: prospective, randomized study of three methods of long-term catheter maintenance. *Crit Care Med.* 1990; 18:1073-9.
34. **Cobb DK, High KP, Sawyer RG, Sable CA, Adams RB, Lindley DA, et al.** A controlled trial of scheduled replacement of central venous and pulmonaryartery catheters. *N Engl J Med.* 1992; 327:1062-8.
35. **Snedecor GW, Cochran WG.** Statistical Methods. Ames, IA: Iowa State Univ Pr; 1978:593.
36. **Maki DG.** Nosocomial bacteremia. An epidemiologic overview. *Am J Med.* 1981; 70:719-32.
37. **Collignon PJ, Munro R, Sorrell TC.** Systemic sepsis and intravenous devices. A prospective survey. *Med J Aust.* 1984; 141:345-8.
38. **Velez LA, Maki DG.** Analysis of risk factors for mortality in nosocomial bloodstream infection. A case-control study [Abstract]. In: Program and Abstracts of the 32nd Interscience Conference on Antimicrobial Agents and Chemotherapy, October 1992, Anaheim, CA. Washington, DC: American Society for Microbiology; 1992:1164.
39. **Gatell JM, Trilla A, Latorre X, Almela M, Mensa J, Moreno A, et al.** Nosocomial bacteremia in a large Spanish teaching hospital: analysis of factors influencing

prognosis. *Rev Infect Dis*. 1988; 10:203-10.

40. Pittet D, Tarara D, Wenzel RP. Nosocomial bloodstream infection in critically ill patients. Excess length of stay, extra costs, and attributable mortality. *JAMA*. 1994; 271:1598-601.

41. Pittet D, Wenzel RP. Nosocomial bloodstream infections in the critically ill. [Letter]. *JAMA*. 1995; 91:1819-20.

42. Gristina AG. Biomaterial-centered infection: microbial adhesion versus tissue integration. *Science*. 1987; 237:1588-95.

43. Vaudaux P, Pittet D, Haeberli A, Huggler E, Nydegger UE, Lew DP, et al. Host factors selectively increase staphylococcal adherence on inserted catheters: a role for fibronectin and fibrinogen or fibrin. *J Infect Dis*. 1989; 160:865-75.

44. Trooskin SZ, Donetz AP, Harvey RA, Greco RS. Prevention of catheter sepsis by antibiotic bonding. *Surgery*. 1985; 97:547-51.

45. Sheretz RJ, Carruth WA, Hampton AA, Byron MP, Solomon DD. Efficacy of antibiotic-coated catheters in preventing subcutaneous *Staphylococcus aureus* infection in rabbits. *J Infect Dis*. 1993; 167:98-106.

46. Gardner JF, Gray KG. Chlorhexidine. In: Block SS, ed. *Disinfection, Sterilization and Preservation*. Philadelphia: Lea & Febiger; 1983:251-70.

47. Freney J, Husson MO, Gavini F, Madier S, Martra A, Izard D, et al. Susceptibilities to antibiotics and antiseptics of new species of the family Enterobacteriaceae. *Antimicrob Agents Chemother*. 1988; 32:873-6.

48. Kahan A, Philippon A, Paul G, Weber S, Richard C, Hazebroucq G, et al. Nosocomia infections by chlorhexidine solution contaminated with *Pseudomonas pickettii* (Biovar VA-I). *J Infect*. 1983; 7:256-63.

49. Stickler DJ, Thomas B, Clayton CL, Chawla JC. Studies of the genetic basis of chlorhexidine resistance. *Br J Clin Pract*. 1983; 25:23-30.

50. Marrie TJ, Costerton JW. Prolonged survival of *Serratia marcescens* in chlorhexidine. *Appl Environ Microbiol*. 1981; 42:1093-102.

51. Aly R, Maibach HI. Effect of antimicrobial soap containing chlorhexidine on the microbial flora of skin. *Appl Environ Microbiol*. 1976; 31:931-5.

52. Okano M, Nomura M, Hata S, Okada N, Sato K, Kitano Y, et al. Anaphylactic symptoms due to chlorhexidine gluconate. *Arch Dermatol*. 1989; 125:50-2.

53. Reynolds NJ, Haman RR. Allergic contact dermatitis from chlorhexidine diacetate in a skin swab. *Contact Dermatitis*. 1990; 22:103-28.

54. Case DE, McAinsh J, Rushton A, Winrow MJ. Chlorhexidine: attempts to detect percutaneous absorption in man. *Chemotherapy*. 1976; 3:367-74.

55. Cowen J, Ellis SH, McAinsh J. Absorption of chlorhexidine from the intact skin of newborn infants. *Arch Dis Child*. 1979; 54:379-83.

56. O'Neill J, Hosmer M, Challop R, Driscoll J, Speck W, Sprunt K. Percutaneous absorption potential of chlorhexidine in neonates. *Current Therapeutic Research*. 1982; 31:485-9.

57. O'Brien CA, Blumer JL, Speck WT, Carr H. Effect of bathing with a 4 percent

chlorhexidine gluconate solution on neonatal bacterial colonization [Abstract]. *J Hosp Infect.* 1984; 5(Suppl A):141.

58. Aoyama H, Yokoo K, Fujii K. Systemic absorption of sulphadiazine, silver sulphadiazine and sodium sulphadiazine through human burn wounds. *Burns.* 1990; 16:163-5.

59. Kulik MI, Wong R, Okarma TB, Falces E, Berkowitz RL. Prospective study of side effects associated with the use of silver sulfadiazine in severely burned patients. *Ann Plast Surg.* 1985; 14:407-19.

60. Degreef H, Dooms-Goossens A. Patch testing with silver sulfadiazine cream. *Contact Dermatitis.* 1985; 12:33-7.

61. Valsecchi R, Foiadelli L, Cainelli T. Is silver sulfadiazine a sensitizer? *Contact Dermatitis.* 1986; 15:45-6.

62. Fuller FW, Engler PE. Leukopenia in non-septic burn patients receiving topical 1% silver sulfadiazine cream therapy: a survey. *J Burn Care Rehabil.* 1988; 9:606-9.

63. Kamal GD, Pfaller MA, Rempe LE, Jebson PJ. Reduced intravascular catheter infection by antibiotic bonding. A prospective, randomized, controlled trial. *JAMA.* 1991; 265:2364-8.

64. Raad I, Darouiche R. Central venous catheters coated with minocycline and rifampin for the prevention of catheter-related bacteremia [Abstract]. In: Program and Abstracts of the 35th Interscience Conference on Antimicrobial Agents and Chemotherapy, October 1995, San Francisco, CA. Washington, DC: American Society for Microbiology; 1995. [Published as Raad I, Darouiche R, Dupuis J, Abi-Said D, Gabrielli A, Hachem R, et al. Central venous catheters coated with minocycline and rifampin for the prevention of catheter-related colonization and bloodstream infections. A randomized, double-blind trial. *Ann Intern Med.* 1997; 127:267-74.].

65. Quesnel LB, Al-Najjar AR, Buddhavudhikrai P. Synergism between chlorhexidine and sulphadiazine. *J Appl Bacteriol.* 1978; 45:397-403.

66. Raad II, Hohn DC, Gilbreath BJ, Suleiman N, Hill LA, Brusio PA, et al. Prevention of central venous catheter-related infections by using maximal sterile barrier precautions during insertion. *Infect Control Hosp Epidemiol.* 1994; 15:231-6.

67. Pearson ML. Guidelines for prevention of intravascular-device-related infections. Hospital Infection Control Practices Advisory Committee. *Infect Control Hosp Epidemiol.* 1996; 17:438-73.